




No.	Plan	Article	Qty
1	P00083261DP	29011	1
2	P00083252DP	29010	1
3	P00083255DP	29009	1
4	XP01795	2993539	2
5	YP06216	27848	1
6	YP06217	27953	1
7	1001680	1001680	2
8	1006999	1006999	2
9	2993761	2993761	2
10	2993786 M5x40	2993786	2
11	1992804	1992804	2
12	1992510	1992510	1
13	1003199	1003199	12
14	P00102003	2998639	4**/0
15	P00095081	2998588	2
16	2998601 SterKnop M8x15	2998601	2

17	1993236	1993236	2
18	X002237	1001688	2
19	2993705BoltM8 x40	2993705	2
20	1001682	1001682	4
21	1992349M8x20	1992349	2
22	P00097039	2998584	4
23	1905365	1998885	2+2**
24	1000972	1000972	4
25	1007000	1007000	8
26	XP01899	GP127	4
27	YP08882	26293	2
28	2993606 Spanring DIN471 15,7-1mm	2993606	2
29	XP02992	GP224	2
30	LA27	1903979	1+1**
31	P00099644	1909002	2 / 0**
32	1912902	1912902	1+1**
33	Label "4" for cable	1995838	1

34	P00045881Bent	2997529	1
35	P00112381	2911819	4**/0
36	XP01906	2993600	2
37	XP01912	2993607	2
38	2993636	2993636	2
39	2993637	2993637	4
40	1006072ScrewP arker4.8x13	1006072	2
41	1001550Washer DIN125M5	1001550	2
42	X004563	1904874	2
43	P00117177	2918919	1
44	P00117221	2919019	1
45	1906684	1906684	6
46	1915944 M8x20	1915944	2

* Not depicted / nie przedstawiono



** 2F option - without legrest motor/ Wersja 2F - bez motoru podnozka

Use only with metal side rails (EP00256 & EP00257),

Do not use with wooden bars (EP00280) and SSR (EP00493)

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS
PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN
WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE,
EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION
FROM VERMEIREN GROUP.

FOOT FRAME ASSEMBLY IZYCO				2991132
Scale: 1:8				Plan No.
Date	2025-01-31	2026-06-16	2026-06-16	E00070499/K.1
Initials	AK	KK	AZ	



Sheet 3/4

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS
PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN
WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE,
EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION
FROM VERMEIREN GROUP.

2991132


Plan No.

E00070499/K.1

FOOT FRAME ASSEMBLY IZYCO

Scale: 1:8	Drafted	Verified	Confirmed
Date	2025-01-31	2026-06-16	2026-06-16
Initials	AK	KK	AZ





2F - use this hole for rastromat mointing

3F - use this hole for rastromat mointing

Sheet 4/4

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS
PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN
WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE,
EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION
FROM VERMEIREN GROUP.

2991132

FOOT FRAME ASSEMBLY IZYCO

Plan No.

E00070499/K.1

Scale:	1:8	Drafted	Verified	Confirmed
Date		2025-01-31	2026-06-16	2026-06-16
Initials		AK	KK	AZ

