




No.	Plan	Article	Qty	20	M6X16 ISO7380	1008459	8	40**	WASHER M4 DIN125 A 4.3X9X0.8	1003303	2	
1	P00084319_DP	29006	1	21#	SCREW M6X12	1007000	8					
2	P00083259_DP	29007	1	22#	XP01899	GP127	4		41*	Label "2" for cable	1995836	1
3	LA27	1901113	1	23#	YP08882	26293	2		42*	Motor cable	1912902	1
4#	XP01795	2993539	2	24#	Spanring DIN471 15,7- 1mm	2993606	2	43**	2998736	2998736	1	
5	YP05894	27847	2									
6	YP06216	27848	1	25#	XP02992	GP224	2					
7	YP06217	27953	1	26	SCREW M5x40	2993786	2					
8	Nut M6 Selflocking	1001680	2	27	Nut M5 DIN986	1992804	2					
9	M6x25 ISO7380	1006998	2	28	WASHER PA6 DIN125A M5	1000972	4					
10	---	2993761	2									
11	Pin D10x50	1998885	2	29	X004563	1904874	2					
12	X001842	G0534	2	30	Distance cap	1001637	2					
13	Tube cap 20x20	1001612	6	31**	P00087487	2998600	2					
14	YP06448	27298	1	32**	M4X10 DIN7991	1907582	2					
15	Pin D8x40	1992510	2	33	XP01906	2993600	2					
16	P00095081	2998588	2	34	XP01912	2993607	2					
17	P00097039	2998584	4	35	Sticker 1	2993636	4					
18	M8X45 ISO7380	1994473	2	36	Sticker 2	2993637	1					
19	NUT M8 SELFLOCKING	1001682	2	37	Washer M5 DIN125	1001550	2					
				38**	Velcro	2992659	2					
				39	Screw Parker 4.8x13	1006072	2					

* not depicted / nie przedstawiono


** Use only on IZYCO

- Use only with metal side rails (EP00256 & EP00257),

Do not use with wooden bars (EP00280) and SSR (EP00493).

Sheet 2/3

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE, EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION FROM VERMEIREN GROUP.				2990966
HEAD SIDE WITH LATHS / IZYCO				Plan No.
				E00073286/I.1
Scale: 1:8	Drafted	Verified	Confirmed	
Date	30.01.2025	16.06.2026	16.06.2026	
Initials	KK	KK	AZ	



NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS
PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN
WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE,
EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION
FROM VERMEIREN GROUP.

2990966

Plan No.

E00073286/1.1

HEAD SIDE WITH LATHS / IZYCO

Scale:	1:8	Drafted	Verified	Confirmed
Date	30.01.2025	16.06.2026	16.06.2026	
Initials	KK	KK	AZ	

