

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE, EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION FROM VERMEIREN GROUP.

### Differential for 2219 with Free Wheel

Scale:	1:5	Drafted	Verified	Confirmed
Date	02/10/2014	2.12.2025	2.12.2025	
Initials	AZ	PM	AZ	

---

Plan No.  
EP00187/I

Material			
No.	Plan	Article	Quantity
1	YP04707	2886600	2
2	XP02365	GP17200	1
3	X002282	1992854	6
4	6x6x50 DIN6885A	GP19000	2
5	XP02371	GP17600	1
6	XP02370	GP17500	1
7	XP01189	2886004	2
8	XP01195	2992221	2
9	XP01197	2992222	2
10	XP01202	2992224	1
11	XP03236	2631200	1
12	XP01283	2992321	1
13	XP01319 (23T)	GP06200	1
	XP03162 (22T)	GP25300	1
14	D10,5 DIN 125	1001558	3
15	M10 DIN 934	1001672	4
16	D15 DIN 6799	1003049	4
17	(Free Wheel -23T)	1008502	1
	XP03160 (Free Wheel -22T)	2996728	1
18	Front Wheel 24"	2992350	1
19	EF00105 Wheel 24" Drive	2994425	2
20**	---	---	-
	Vintage	1907965	3
21*	---	1902733	3
22	M8x25	1990501	2
23	M5x16	1006112	6
24	---	1994568	6
25	D8,4 DIN 9021	1996574	2
26	M10 DIN 493B	1903216	1
27	M6x90 DIN 912	1911103 (01,02,03)	4
	M6x110 DIN 912	1911104 (06)	4
	M6x100 DIN 912	1998869 (07)	4
28	M6 DIN 127B	1004414	4


\* not depicted

\*\* option

Material			
No.	Plan	Article	Quantity
29	Chain	2993196	1,4 [m]
30	---	G040600	1
31	XP01259	2992306 (01)	1
	XP01671	2993508 (02)	1
	XP01672	2993509 (03)	1
	XP01673	2993510 (06)	1
	XP01674	2993511 (07)	1

2/2

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE, EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION FROM VERMEIREN GROUP.

Differential for 2219 with Free Wheel				---
				Plan No. EP00187/I
Scale: 1:5	Drafted	Verified	Confirmed	
Date	02/10/2014	2.12.2025	2.12.2025	
Initials	AZ	PM	AZ	