

No.	Plan	Article	Quantity
1	---	L000068	0,5
2	---	1910488	1
3	---	1994180	2
4	---	1901731	2
5	---	1007000	2
6	---	1004414	2
7	---	1001680	2
8	---	L000007	2
9	---	L000001	4
10	---	T001009	0,1


Turios Battery Harness Set 2 is used for:  
 Turios with lights and without electrical seat and backrest angle adjustment

TURIOS BATTERY HARNESS SET 2 Sheet 1/2

Scale:	1:1	Drafted	Verified	Confirmed
Date	04/11/2021	04/11/2021	04/11/2021	04/11/2021
Initials	KS	MVG	TK	

---

Plan No.  
E00008120/A



12

3

9

4

12

14

15

13

11

330mm

550mm

800mm

No.	Plan	Article	Quantity
11	---	1994573	1
12	---	1901755	2
13	---	H005006	0,33
14	---	1907749	0,8
15	---	2990811	2

TURIOS BATTERY HARNESS SET 2 Sheet2/2

Scale:	1:1	Drafted	Verified	Confirmed
Date	04/11/2021	04/11/2021	04/11/2021	04/11/2021
Initials	KS	MVG	TK	

---
Plan No.
E00008120/A
